



# GeloPure Q

**Ion Exchange Chromatography Media**

## **Technical Datasheet**

## GeloPure Q

GeloPure Q strong anion exchanger is a surface-modified, highly cross-linked medium with dextran scaffold structure. It delivers high dynamic binding capabilities and maintains stability even at elevated flow rates. The optimized media's excellent performance presents a substantial chance to boost the throughput of downstream purification. In this ion exchange resin, the Quaternary amine functional group is attached to the Cross-linked agarose beads, which remains charged and maintains consistently high capacity over the entire working range, pH 2 to 12.




### GeloPure Q base resin

Cross-linked spherical agarose beads are prepared by unique process from natural polysaccharides. Primarily composed of agarobiose repeating units (alternating D-galactose and 3,6-anhydro-L-galactose), finally beads are highly hydrophilic, minimize the non-specific binding. GeloPure Q, made from cross-linked agarose 4% beads, benefits from this structure by offering efficient flow velocity for large biomolecules. This helps improve performance when purifying large proteins and other biological substances.

### Pressure-flow and Characteristics of GeloPure Q

GeloPure Q media consist of highly cross-linked agarose beads. It gives high binding capacities even in high flowrates. It allows flow rates of 600 cm/h through at 1 bar (14.5 psi, 0.1 MPa), enabling quick separation steps—especially useful in the early stages of purification

when fast processing is important. During washing and equilibration, the flow rate can be increased up to 750 cm/h for even faster operation. The basic characteristics of GeloPure Q weak anion exchange chromatography media are shown in Table 1.



Ion Exchange Type		Strong anion exchange
Ligand name		-N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>
Base Matrix		Highly Cross-linked agarose 4%
Particle Size (µm)		40 µm to 165 µm
pH Working Range		2 to 12
Operating Pressure		Up to 2 Bar (0.2 Mpa)
Ion Exchange Capacity (meq / ml-gel)		0.10-0.24
Chemical Stability		0.5N NaOH
Dynamic Binding Capacity (mg/ml)	HSA	>100
	human-γ-globulin	35
Supplied		Suspension in 20% Ethanol

**Table 1.** Basic characteristics of GeloPure Q strong anion exchange chromatography media

## Dynamic Binding Capacities of GeloPure Q

GeloPure Q have high efficiency in mass transfer and excellent Dynamic Binding Capacities, particularly for HSA (Human Normal Albumin) and large biomolecules like Immunoglobulins (IgG).

Because of these special qualities, GeloPure Q media can be used in downstream processes in the purification of biopharmaceuticals. With good reliability, protein separation has been scaled up using GeloPure Q from laboratory size and no discernible variations were seen in the purity of single peak at the time of protein separation.

Column: HiScale 16/20, 5 cm bed height, 10 ml bed volume

Flow rate: 1ml/min

Sample: 20mg/ml

Start Buffer: 30mM Phosphate buffer (pH 7.8-8.2) for HNA

50mM Tris-HCL (pH 9.5) for IgG

Elution Buffer: 30mM Sodium Acetate (pH 4.6) for HNA

50mM Tris-HCL+1M NaCl (pH 9.5) for IgG

## Chemical Stability and Cleaning-In-Place

GeloPure Q, a cross-linked spherical GeloPure beads, offer strong chemical and physical stability, enabling effective clean-in-place (CIP) and sanitization protocols. This ensures high protein recovery over multiple cycles, prevents microbial growth, and supports hygienic, cost-effective purification which is the key considerations for preparative applications.

Cleaning-In-Place of GeloPure Q done with 0.5N NaOH and sometimes with Orthophosphoric acid. Most of the contaminated material should be removed from CIP with routine washing with 5 CV of 0.25N to 0.5N sodium hydroxide; however, extremely hydrophobic molecules may bind so firmly that they must be eluted using powerful detergents or organic solvents, such as 70% ethanol or 30% isopropyl alcohol.

After use of chromatographic media, it should be stored in 20% ethanol at room temperature.